

UW



UX



BW-K



BX-K



UW Series



BW-K plus  
Medium-size Animal Bucket



BW-K plus  
Large Animal Bucket

Electronic Balances for Weighing Animals

**Dedicated software brings quick and reliable results of live animal weighing.**

Upon removing the weighed animal, the weight of remaining feces and output are automatically tared off, and the balance presents zero for the next animal.

Order a standard balance plus one of the applicable animal bucket sets.

Model	Catalog No.	Balance Series	Reduced Capacity (approx.)
Small Animal Bucket set	321-62150	UW/UX (Capacity 2200g or more)	Bottom 110dia, Top 200dia, Height 130
Medium-size Animal Bucket set <sup>*1</sup>	321-62545-03	BX-K	Bottom 305 x 215, Top 377 x 245, Height 215
	321-62545-04	BW-K	Bottom 305 x 215, Top 377 x 245,

			Height 215
Large Animal Bucket set <sup>*2</sup>	321-62545- 01	BX-K (Capacity 22kg or more)	Bottom 335 x 245, Top 445 x 395, Height 345
	321-62545- 02	BW-K (Capacity 22kg or more)	Bottom 335 x 245, Top 445 x 395, Height 345

\*1 Capacity is reduced about 2kg

\*2 Capacity is reduced about 6kg

For Research Use Only. Not for use in diagnostic procedures.